

Chlorine Disinfection



Artikelnr.: 60152



Artikelnr.:60132

Broad-spectrum disinfection
Safe-to-use
Always correct dosing
Long-term storage stability

Wet Wipe Chlorine Disinfection Wipes are impregnated with a granulated chlorine product
By using Wet Wipe Chlorine Disinfection you will experience a safe work procedure and achieve an optimal disinfection. A pack of products consists of 10 wipes impregnated with a stable chlorine granulate and a plastic bag filled with water.

Instruction for use
The product is activated by hitting the pack, which cracks the bag, such that water burst into the wipes. When the granulate dissolves, it is turned into a solution of hypochlorous acid with an active chlorine content of 2000 ppm. Wipes are activated and ready to use after 7 minutes. Can be used for disinfection of hard surfaces that can withstand chlorine. Wipe of surfaces with a clean damp cloth after disinfection. Use carefully on lacquered surfaces as discoloration may occur. Do not use for disinfection of natural rubber materials.

EfficacyEfficacy
Effective against vegetative bacteria, mycobacteria, yeast, mould fungus, viruses and certain bacterial spores. Recommended for terminal disinfection and for disinfection of surroundings of patients in isolation with *Clostridium difficile* and viruses.

Stability
Products has 2 years of stability before activation. Activated ready to use wipes retain their disinfection efficacy up to 7 days after activation.

The wipes are produced in two sizes

Item no.: 60152: Chlorine Disinfection Maxi
Karton á 24 pakker
1 pakke indeholder 5 klude, 43x30 cm.

Item no.: 60132: Chlorine Disinfection Mini
Carton of 20 packets
1 packet contains 10 wipes, 20x30 cm

Content
Chlorine granules: Sodium dichloroisocyanurate dihydrate 100%.
Activated ready to use wipes: Active chlorine > 2000 ppm.

Always read the product labeling before use.
Safety datasheet for the activated ready to use product is available.

Produced in Germany.

Wet Wipe A/S
Tel: +45 70 266 244
info.wetwipe@schuelke.com
www.wetwipe.eu

CVR/VAT no: DK 30207726

DK/EN, Ver4
Maj 2024

wetwipe®
by schülke →

The disinfection efficacy have been tested under clean conditions with the following results.

Minimum application time for pass of EN test

Microorganism	EN test	Active chlorine	Time
Bacteria: <i>Staphylococcus aureus</i> , <i>Enterococcus hirae</i> , <i>Escherichia coli</i> , <i>Pseudomonas aeruginosa</i>	EN 13727 EN 16615	1000 ppm 1200 ppm	1 min 5 min
Yeast fungus: <i>Candida albicans</i>	EN 13624 EN 16615	1000 ppm 1200 ppm	1 min 5 min
Mould fungus spores: <i>Aspergillus brasiliensis</i>	EN 13624	2000 ppm	15 min
Virus: <i>Norovirus</i> <i>Adenovirus</i> , <i>poliovirus</i>	EN 14476	500 ppm	5 min
Tuberculosis: <i>Mycobacterium terrae</i> , <i>Mycobacterium avium</i>	EN 14348	2000 ppm	5 min
Bacterial spores: <i>Clostridium difficile</i>	EN 17126	1000 ppm	5 min